



## Sessional II(odd semester) Examination, October 2025

Roll no.....

Name of the Course and semester: B. Pharm. Isem

Name of the Paper: Pharmaceutical Analysis-I

Paper Code: BP-102T

Time: 1.5-hour

Maximum Marks: 30

**Note:**

- (i) This question paper contains three sections.
- (ii) All the questions are compulsory.

### Section-A

**Q1. Multiple choice questions (attempt all questions)**

10×1= 10 Marks

1. EDTA titrations mainly belong to which class? CO3
  - a. Acid-base
  - b. Redox
  - c. Complexometric
  - d. Precipitation
  
2. A masking agent is used to CO3
  - a) Precipitate the analyte
  - b) Form a stable, non-titrable complex with an interfering ion
  - c) Oxidize the analyte
  - d) Change the oxidation state of EDTA
  
3. Demasking refers to: CO3
  - a) Back-titration step
  - b) Releasing a previously masked ion to allow its subsequent titration
  - c) Changing indicator midway
  - d) Increasing ionic strength
  
4. Gravimetric analysis is based on: CO3
  - a) Volume
  - b) Potential
  - c) Mass of a compound of known stoichiometry
  - d) Color intensity
  
5. Co-precipitation refers to: CO3
  - a) Precipitation of analyte only
  - b) Carry-down of impurities with precipitate by inclusion/occlusion/adsorption
  - c) Loss of analyte during washing
  - d) Volatilization upon ignition

6. In diazotisation, the reagent is typically: CO3
- Sodium nitrite in acid at 0–5 °C
  - Potassium permanganate at 60 °C
  - Sodium thiosulfate in base
  - Ferric chloride in ethanol
7. A common endpoint detection in diazotisation titration is: CO3
- Internal azo dye formation always
  - External starch-iodide paper (blue color in excess nitrite)
  - Ferric alum red color
  - Chromate red-brown
8. To minimize co-precipitation, one should: CO 3
- Precipitate hot, from dilute solution, with slow addition and digestion
  - Precipitate cold and fast
  - Use highly supersaturated conditions
  - Avoid washing
9. Post-precipitation means: CO 3
- Formation of a new precipitate after the main precipitate has formed.
  - Re-dissolution of precipitate
  - Breakage of crystals by stirring
  - Adsorption of indicator
10. The correct sequence in precipitation gravimetry is: CO3
- Filtration → precipitation → digestion → ignition → weighing
  - Precipitation → digestion (aging) → filtration/wash → drying/ignition → weighing
  - Digestion → ignition → precipitation → filtration
  - Precipitation → ignition → wash → weighing

### Section B

- Q2. Short type (attempt any two questions)** 2×5= 10Marks
- Difference between masking and demasking agents. CO 3
  - Classify metal ion indicators. CO 3
  - Explain the process of gravimetry. CO 3

### Section-C

- Q3. Long type (attempt any one questions)** 1×10= 10 Marks
- Classify different types of precipitation titrations and their applications. CO 3
  - Discuss Complexometric titrations in detail and its application in estimation of magnesium sulphate. CO 3